

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MIA Polyclonal Antibody

多克隆抗体

产品货号: PA04957

第三版

**描述:**MIA (Melanoma Inhibitory Activity) is a Protein Coding gene. Diseases associated with MIA include Melanoma and Skin Melanoma. Among its related pathways are Neural Crest Differentiation. GO annotations related to this gene include growth factor activity. An important paralog of this gene is MIA-RAB4B.

**抗原:**Synthetic peptide from human protein

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**All malignant melanoma cell lines tested and infrequently in glioma cell lines.