

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MGST2 Polyclonal Antibody

多克隆抗体

产品货号: PA04963

第三版

描述:The MAPEG (Membrane Associated Proteins in Eicosanoid and Glutathione metabolism) family consists of six human proteins, several of which are involved in the production of leukotrienes and prostaglandin E, important mediators of inflammation. This gene encodes a protein which catalyzes the conjugation of leukotriene A4 and reduced glutathione to produce leukotriene C4. Alternatively spliced transcript variants encoding different isoforms have been identified in this gene. MGST2 (Microsomal Glutathione S-Transferase 2) is a Protein Coding gene. Among its related pathways are Platinum drug resistance and Glutathione metabolism. GO annotations related to this gene include enzyme activator activity and glutathione peroxidase activity. An important paralog of this gene is LTC4S.

抗原:Synthesized peptide derived from the Internal region of human MGST2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Liver, spleen, skeletal muscle, heart, adrenals, pancreas, prostate, testis, fetal liver, and fetal spleen. Very low expression in lung, brain, placenta and bone marrow.