

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD52 Polyclonal Antibody

多克隆抗体

产品货号: PA04977

第三版

描述:CD52 (CD52 Molecule) is a Protein Coding gene. Diseases associated with CD52 include Relapsing-Remitting Multiple Sclerosis and Chronic Lymphocytic Leukemia. Among its related pathways are Post-translational modification- synthesis of GPI-anchored proteins and Metabolism of proteins.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: