

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat/Monkey SPON2 Polyclonal Antibody

多克隆抗体

产品货号: PA04978

第三版

描述:SPON2 (Spondin 2) is a Protein Coding gene. Among its related pathways are Phospholipase-C Pathway and O-glycosylation of TSR domain-containing proteins. GO annotations related to this gene include antigen binding and lipopolysaccharide binding. An important paralog of this gene is SPON1.

抗原:Synthesized peptide derived from the Internal region of human Mindin

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in normal lung tissue but not in lung carcinoma cell lines.