

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human CD300A Polyclonal Antibody

多克隆抗体

产品货号: PA04985

第三版

**描述:**CD300A (CD300a Molecule) is a Protein Coding gene. Diseases associated with CD300A include Kabuki Syndrome 1. Among its related pathways are Immunoregulatory interactions between a Lymphoid and a non-Lymphoid cell and Innate Immune System. GO annotations related to this gene include phosphatidylserine binding and signaling receptor activity. An important paralog of this gene is CD300C. This gene encodes a member of the CD300 glycoprotein family of cell surface proteins found on leukocytes involved in immune response signaling pathways. This gene is located on chromosome 17 in a cluster with all but one of the other family members. Multiple transcript variants encoding different isoforms have been found for this gene.

**抗原:**Synthetic peptide from human protein

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed not only by natural killer (NK) cells but also by T-cell subsets, B-cells, dendritic cells, mast cells, granulocytes and monocytes.