

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MKP-1 Polyclonal Antibody

多克隆抗体

产品货号: PA04991

第三版

描述: Dual specificity phosphatase that dephosphorylates MAP kinase ERK2 on both 'Thr-183' and 'Tyr-185'.

抗原: Synthesized peptide derived from human MKP-1 around the non-phosphorylation site of Ser359.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed at high levels in the lung, liver placenta and pancreas. Moderate levels seen in the heart and skeletal muscle. Lower levels found in the brain and kidney.