

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MKP1/2 Polyclonal Antibody

多克隆抗体

产品货号: PA04995

第三版

描述:The expression of DUSP1 gene is induced in human skin fibroblasts by oxidative/heat stress and growth factors. It specifies a protein with structural features similar to members of the non-receptor-type protein-tyrosine phosphatase family, and which has significant amino-acid sequence similarity to a Tyr/Ser-protein phosphatase encoded by the late gene H1 of vaccinia virus. The bacterially expressed and purified DUSP1 (Dual Specificity Phosphatase 1) is a Protein Coding gene. Diseases associated with DUSP1 include Vaccinia and Multiple Cranial Nerve Palsy. Among its related pathways are RET signaling and Serotonergic synapse. GO annotations related to this gene include phosphatase activity and protein tyrosine/serine/threonine phosphatase activity. An important paralog of this gene is DUSP4. DUSP4 (Dual Specificity Phosphatase 4) is a Protein Coding gene. Diseases associated with DUSP4 include Echolalia and Chronic Laryngitis. Among its related pathways are RET signaling and Activated TLR4 signalling. GO annotations related to this gene include phosphatase activity and protein tyrosine/serine/threonine phosphatase activity. An important paralog of this gene is DUSP1.

抗原:Synthesized peptide derived from human MKP-1/2 around the non-phosphorylation site of S296/318.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed at high levels in the lung, liver placenta and pancreas. Moderate levels seen in the heart and skeletal muscle. Lower levels found in the brain and kidney.1 Publication