

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Granzyme M Polyclonal Antibody

多克隆抗体

产品货号: PA05007

第三版

描述:GZMM (Granzyme M) is a Protein Coding gene. Diseases associated with GZMM include Adrenal Insufficiency, Congenital, With 46Xy Sex Reversal, Partial Or Complete and Congenital Adrenal Insufficiency. Among its related pathways are Innate Immune System and Granzyme Pathway. GO annotations related to this gene include serine-type endopeptidase activity and endopeptidase activity. An important paralog of this gene is CFD. Human natural killer (NK) cells and activated lymphocytes express and store a distinct subset of neutral serine proteases together with proteoglycans and other immune effector molecules in large cytoplasmic granules. These serine proteases are collectively termed granzymes and include 4 distinct gene products: granzyme A, granzyme B, granzyme H, and the protein encoded by this gene, granzyme M. Two transcript variants encoding different isoforms have been found for this gene.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μg/ml.

特异性:Highly and constitutively expressed in activated natural killer (NK) cells.