

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse EYA1/4 Polyclonal Antibody

多克隆抗体

产品货号: PA05008

第三版

描述:EYA1 (EYA Transcriptional Coactivator And Phosphatase 1) is a Protein Coding gene. Diseases associated with EYA1 include Otofaciocervical Syndrome and Branchiootoc Syndrome 1. Among its related pathways are DNA Double-Strand Break Repair and Transcriptional misregulation in cancer. GO annotations related to this gene include RNA binding and protein tyrosine phosphatase activity. An important paralog of this gene is EYA4. EYA4 (EYA Transcriptional Coactivator And Phosphatase 4) is a Protein Coding gene. Diseases associated with EYA4 include Cardiomyopathy, Dilated, 1J and Deafness, Autosomal Dominant 10. Among its related pathways are DNA Double-Strand Break Repair and DNA Double Strand Break Response. GO annotations related to this gene include protein tyrosine phosphatase activity. An important paralog of this gene is EYA1.

抗原:Synthetic peptide from human protein

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: