
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MMP1 Polyclonal Antibody

多克隆抗体

产品货号: PA05016

第三版

描述: Cleaves collagens of types I, II, and III at one site in the helical domain. Also cleaves collagens of types VII and X. In case of HIV infection, interacts and cleaves the secreted viral Tat protein, leading to a decrease in neuronal Tat's mediated neurotoxicity.

抗原: Synthesized peptide derived from the C-terminal region of human MMP-1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: