

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MN1 Polyclonal Antibody

多克隆抗体

产品货号: PA05024

第三版

**描述:** Meningioma 1 (MN1) contains two sets of CAG repeats. It is disrupted by a balanced translocation (4;22) in a meningioma, and its inactivation may contribute to meningioma pathogenesis. MN1 (MN1 Proto-Oncogene, Transcriptional Regulator) is a Protein Coding gene. Diseases associated with MN1 include Meningioma, Familial and Mn1.

**抗原:** Synthesized peptide derived from the Internal region of human MN1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:**