
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MKNK2 Polyclonal Antibody

多克隆抗体

产品货号: PA05031

第三版

描述:MKNK2 (MAP Kinase Interacting Serine/Threonine Kinase 2) is a Protein Coding gene. Among its related pathways are Type I Interferon Signaling Pathways and TGF-Beta Pathway. GO annotations related to this gene include transferase activity, transferring phosphorus-containing groups and protein tyrosine kinase activity. An important paralog of this gene is MKNK1.

抗原:Synthesized peptide derived from the Internal region of human Mnk2.

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Ubiquitously expressed in all tissues examined. Isoform 2 is expressed at higher levels in the ovary than is isoform 1.