



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NUTM1 Polyclonal Antibody

多克隆抗体

产品货号: PA05038

第三版

描述:NUTM1 (NUT Midline Carcinoma Family Member 1) is a Protein Coding gene. Diseases associated with NUTM1 include Nut Midline Carcinoma and Extragenital Nonseminomatous Germ Cell Tumor. Among its related pathways are Chromatin Regulation / Acetylation. An important paralog of this gene is NUTM2D.

抗原:Synthesized peptide derived from human NUT

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Specifically expressed in testis.