

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MOV10L1 Polyclonal Antibody

多克隆抗体

产品货号: PA05044

第三版

**描述:** This gene is similar to a mouse gene that encodes a putative RNA helicase and shows testis-specific expression. Multiple transcript variants encoding different isoforms have been found for this gene. MOV10L1 (Mov10 RISC Complex RNA Helicase Like 1) is a Protein Coding gene. Among its related pathways are Gene Expression and Mitotic Prophase. GO annotations related to this gene include RNA binding and ATP-dependent RNA helicase activity. An important paralog of this gene is MOV10.

**抗原:** Synthesized peptide derived from the Internal region of human MOV10L1

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5 μ g/ml.

**特异性:** Isoform 1: Specifically expressed in testis.