

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MRGPRX3 Polyclonal Antibody

多克隆抗体

产品货号: PA05060

第三版

描述: This gene encodes a member of the mas-related/sensory neuron specific subfamily of G protein coupled receptors. The encoded protein may be involved in sensory neuron regulation and in the modulation of pain. MRGPRX3 (MAS Related GPR Family Member X3) is a Protein Coding gene. GO annotations related to this gene include G-protein coupled receptor activity. An important paralog of this gene is MRGPRX1.

抗原: Synthesized peptide derived from the Internal region of human MRGX3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Uniquely localized in a subset of small dorsal root and trigeminal sensory neurons.