

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat TNFAIP8L3 Polyclonal Antibody

多克隆抗体

产品货号: PA05065

第三版

描述: TNFAIP8L3 (TNF Alpha Induced Protein 8 Like 3) is a Protein Coding gene. Among its related pathways are Metabolism and Glycerophospholipid biosynthesis. An important paralog of this gene is TNFAIP8.

抗原: Synthesized peptide derived from human TNFAIP8L3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed (at protein level) (PubMed:25479791). Highly expressed in most carcinoma cell lines (PubMed:25479791, PubMed:25242044).