

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ABCC3 Polyclonal Antibody

多克隆抗体

产品货号: PA05066

第三版

描述: May act as an inducible transporter in the biliary and intestinal excretion of organic anions. Acts as an alternative route for the export of bile acids and glucuronides from cholestatic hepatocytes.

抗原: Synthesized peptide derived from the Internal region of human MRP3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Mainly expressed in the liver. Also expressed in small intestine, colon, prostate, testis, brain and at a lower level in the kidney.