

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat DLX2 Polyclonal Antibody

多克隆抗体

产品货号: PA05072

第三版

描述:DLX2 (Distal-Less Homeobox 2) is a Protein Coding gene. Diseases associated with DLX2 include Tooth Agenesis and Bamforth-Lazarus Syndrome. Among its related pathways are Sudden Infant Death Syndrome (SIDS) Susceptibility Pathways and Preimplantation Embryo. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and chromatin binding. An important paralog of this gene is DLX5. Many vertebrate homeo box-containing genes have been identified on the basis of their sequence similarity with Drosophila developmental genes. Members of the Dlx gene family contain a homeobox that is related to that of Distal-less (Dll), a gene expressed in the head and limbs of the developing fruit fly. The Distal-less (Dlx) family of genes comprises at least 6 different members, DLX1-DLX6. The DLX proteins are postulated to play a role in forebrain and craniofacial development. This gene is located in a tail-to-tail configuration with another member of the gene family on the long arm of chromosome 2.

抗原: Synthesized peptide derived from human Dlx-2

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: