

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human COL19A1 Polyclonal Antibody

多克隆抗体

产品货号: PA05082

第三版

描述:COL19A1 (Collagen Type XIX Alpha 1 Chain) is a Protein Coding gene. Among its related pathways are Phospholipase-C Pathway and Degradation of the extracellular matrix. GO annotations related to this gene include extracellular matrix structural constituent and protein binding, bridging. An important paralog of this gene is COL22A1. This gene encodes the alpha chain of type XIX collagen, a member of the FACIT collagen family (fibril-associated collagens with interrupted helices). Although the function of this collagen is not known, other members of this collagen family are found in association with fibril-forming collagens such as type I and II, and serve to maintain the integrity of the extracellular matrix. The transcript produced from this gene has an unusually large 3' UTR which has not been completely sequenced.

抗原:Synthesized peptide derived from human COL19A1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Localized to vascular, neuronal, mesenchymal, and some epithelial basement membrane zones in umbilical cord.