

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human OASL Polyclonal Antibody

多克隆抗体

产品货号: PA05116

第三版

**描述:**OASL encodes an interferon-induced protein that is a member of the 2-5A synthetase family. 2-prime,5-prime oligoadenylates (2-5As) bind to and activate RNase L, resulting in general RNA degradation and consequent inhibition of protein synthesis. 2-5As are produced by a well-conserved family of interferon-induced enzymes, the 2-5A synthetases or OASs

**抗原:**Synthesized peptide derived from human OASL

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in most tissues, with the highest levels in primary blood Leukocytes and other hematopoietic system tissues, colon, stomach and to some extent in testis.