

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Phospho-RIPK2 (Ser176) Polyclonal Antibody

多克隆抗体

产品货号: PA05131

第三版

描述: This gene encodes a member of the receptor-interacting protein (RIP) family of serine/threonine protein kinases. The encoded protein contains a C-terminal caspase activation and recruitment domain (CARD), and is a component of signaling complexes in both the innate and adaptive immune pathways. It is a potent activator of NF-kappaB and inducer of apoptosis in response to various stimuli. RIPK2 (Receptor Interacting Serine/Threonine Kinase 2) is a Protein Coding gene. Diseases associated with RIPK2 include Lymphoproliferative Syndrome, X-Linked, 2 and Human Granulocytic Anaplasmosis. Among its related pathways are IL-1 signaling pathway and Nucleotide-binding domain, leucine rich repeat containing receptor (NLR) signaling pathways. GO annotations related to this gene include protein homodimerization activity and protein kinase activity. An important paralog of this gene is RIPK4.

抗原: Synthesized phospho derived from human RIPK2 (Phospho-Ser176)

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) , 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Detected in heart, brain, placenta, lung, peripheral blood leukocytes, spleen, kidney, testis, prostate, pancreas and lymph node.