

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MSH3 Polyclonal Antibody

多克隆抗体

产品货号: PA05133

第三版

**描述:**The protein encoded by this gene forms a heterodimer with Human, Mouse SH2 to form MutS beta, part of the post-replicative DNA mismatch repair system. MutS beta initiates mismatch repair by binding to a mismatch and then forming a complex with Human, Mouse MSH2 alpha heterodimer. This gene contains a polymorphic 9 bp tandem repeat sequence in the first exon. The repeat is present 6 times in the reference genome sequence and 3-7 repeats have been reported. Defects in this gene are a cause of susceptibility to endometrial cancer.

**抗原:**Synthesized peptide derived from the N-terminal region of human MSH3

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5 μ g/ml.

**特异性:**