

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SIN3B Polyclonal Antibody

多克隆抗体

产品货号: PA05134

第三版

描述: Acts as a transcriptional repressor. Interacts with MXI1 to repress MYC responsive genes and antagonize MYC oncogenic activities. Interacts with MAD-MAX heterodimers by binding to MAD. The heterodimer then represses transcription by tethering SIN3B to DNA. Also forms a complex with FOXK1 which represses transcription (By similarity).

抗原: Synthesized peptide derived from the Internal region of human Sin3B

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: