

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat STK24 Polyclonal Antibody

多克隆抗体

产品货号: PA05138

第三版

描述: This gene encodes a serine/threonine protein kinase that functions upstream of mitogen-activated protein kinase (MAPK) signaling. The encoded protein is cleaved into two chains by caspases; the N-terminal fragment (MST3/N) translocates to the nucleus and promotes programmed cells death. There is a pseudogene for this gene on chromosome X. Alternative splicing results in multiple transcript variants.

抗原: Synthesized peptide derived from the C-terminal region of human MST-3

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Isoform A is ubiquitous. Isoform B is expressed in brain with high expression in hippocampus and cerebral cortex.