

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MUC13 Polyclonal Antibody

多克隆抗体

产品货号: PA05174

第三版

描述: Epithelial mucins, such as MUC13, are a family of secreted and cell surface glycoproteins expressed by ductal and glandular epithelial tissues (

抗原: Synthesized peptide derived from the C-terminal region of human Mucin 13

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Highly expressed in epithelial tissues, particularly those of the gastrointestinal and respiratory tracts, such as large intestine and trachea, followed by kidney, small intestine, appendix and stomach.