

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CK-5 Polyclonal Antibody

多克隆抗体

产品货号: PA05177

第三版

描述:In normal tissues, CK5 is expressed in the basal cells of squamous epithelium and ductal epithelium, and part of the germinal cells of the squamous epithelium, myoepithelial cells and mesothelial cells. While it is almost not expressed in the vast majority of glandular epithelial cells. Therefore, it can be used for differential diagnosis of squamous cell carcinoma and adenocarcinoma, mesothelioma and adenocarcinoma. And it can also be used for differential diagnosis of mammary gland, salivary gland, benign and malignant tumor derived from prostate epithelial.

抗原:Synthesized peptide derived from human Keratin 5 Polyclonal

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: