

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat PLD4 Polyclonal Antibody

多克隆抗体

产品货号: PA05205

---

第三版

**描述:**PLD4 (Phospholipase D Family Member 4) is a Protein Coding gene. Among its related pathways are Role of phospholipids in phagocytosis and Innate Immune System. GO annotations related to this gene include phosphatidylinositol phospholipase C activity and N-acylphosphatidylethanolamine-specific phospholipase D activity. An important paralog of this gene is PLD3.

**抗原:**Synthesized peptide derived from human PLD4 Polyclonal

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**