

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PKMYT1 Polyclonal Antibody

多克隆抗体

产品货号: PA05206

第三版

描述: This gene encodes a member of the serine/threonine protein kinase family. The encoded protein is a membrane-associated kinase that negatively regulates the G2/M transition of the cell cycle by phosphorylating and inactivating cyclin-dependent kinase 1. The activity of the encoded protein is regulated by polo-like kinase 1. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.

抗原: Synthesized peptide derived from human Myt 1 around the non-phosphorylation site of Ser83.

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4)，1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: