

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MZF1 Polyclonal Antibody

多克隆抗体

产品货号: PA05208

第三版

描述:MZF1 (Myeloid Zinc Finger 1) is a Protein Coding gene. Diseases associated with MZF1 include Arrhythmogenic Right Ventricular Dysplasia 1. Among its related pathways are Ectoderm Differentiation. GO annotations related to this gene include transcription factor activity, sequence-specific DNA binding and RNA polymerase II core promoter proximal region sequence-specific DNA binding. An important paralog of this gene is ZSCAN10.

抗原:Synthesized peptide derived from the N-terminal region of human MZF-1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Preferentially expressed in differentiating myeloid cells. Detected in osteoblasts.