
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human N4BP1 Polyclonal Antibody

多克隆抗体

产品货号: PA05212

第三版

描述:N4BP1 (NEDD4 Binding Protein 1) is a Protein Coding gene. GO annotations related to this gene include RNA binding. An important paralog of this gene is KHNYN.

抗原:Synthesized peptide derived from the Internal region of human N4BP1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in heart, lung, brain, liver, skeletal muscle, pancreas, kidney, spleen, testis and ovary.