

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SLC22A6 Polyclonal Antibody

多克隆抗体

产品货号: PA05220

第三版

描述:SLC22A6 (Solute Carrier Family 22 Member 6) is a Protein Coding gene. Diseases associated with SLC22A6 include N-Acetylglutamate Synthase Deficiency. Among its related pathways are Transport of glucose and other sugars, bile salts and organic acids, metal ions and amine compounds and Methotrexate Pathway, Pharmacokinetics. GO annotations related to this gene include protein homodimerization activity and ion transmembrane transporter activity. An important paralog of this gene is SLC22A8. The protein encoded by this gene is involved in the sodium-dependent transport and excretion of organic anions, some of which are potentially toxic. The encoded protein is an integral membrane protein and may be localized to the basolateral membrane. Four transcript variants encoding four different isoforms have been found for this gene.

抗原:Synthesized peptide derived from human OAT1 Polyclonal

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: