
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat AZI2 Polyclonal Antibody

多克隆抗体

产品货号: PA05223

第三版

描述: AZI2, or NAP1, contributes to the activation of NF κ B (see MIM 164011)-dependent gene expression by activating IKK-related kinases, such as NAK. AZI2 (5-Azacytidine Induced 2) is a Protein Coding gene. Among its related pathways are RIG-I/MDA5 mediated induction of IFN-alpha/beta pathways.

抗原: Synthesized peptide derived from the N-terminal region of human Nap1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Widely expressed. Abundant expression seen in the pancreas and testis.