

本公司提供的电版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SLC24A4 Polyclonal Antibody

多克隆抗体

产品货号: PA05231

第三版

描述: This gene encodes a member of the potassium-dependent sodium/calcium exchanger protein family. Alternative splicing results in multiple transcript variants. Transports 1 Ca²⁺ and 1 K⁺ in exchange for 4 Na⁺. Controls the rapid response termination and proper regulation of adaptation in olfactory sensory neurons (OSNs) which subsequently influences how odor information is encoded and perceived.

抗原: Synthesized peptide derived from the Internal region of human NCKX4

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Expressed abundantly in all regions of the brain, aorta, lung and thymus. Expressed at lower levels in the stomach and intestine.