

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MT-ND5 Polyclonal Antibody

多克隆抗体

产品货号: PA05233

第三版

**描述:**Core subunit of the mitochondrial membrane respiratory chain NADH dehydrogenase (Complex I) that is believed to belong to the minimal assembly required for catalysis. Complex I functions in the transfer of electrons from NADH to the respiratory chain. The immediate electron acceptor for the enzyme is believed to be ubiquinone.

**抗原:**Synthesized peptide derived from the Internal region of human ND5

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**