
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat/Monkey NEDD8 Polyclonal Antibody

多克隆抗体

产品货号: PA05242

第三版

描述:NEDD8 (neural precursor cell expressed, developmentally down-regulated 8) is a protein-coding gene. Diseases associated with NEDD8 include glomovenous malformation. GO annotations related to this gene include ubiquitin protein ligase binding.

抗原:Synthesized peptide derived from the N-terminal region of human Neddylin

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Highly expressed in heart, skeletal muscle, spleen, thymus, prostate, testis, ovary, colon and leukocytes.