

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NEIL1 Polyclonal Antibody

多克隆抗体

产品货号: PA05246

第三版

描述: This gene is a member of the Nei endonuclease VIII-like gene family which encodes DNA glycosylases. The encoded enzyme participates in the DNA repair pathway by initiating base excision repair by removing damaged bases, primarily oxidized pyrimidines. Multiple transcript variants encoding different isoforms have been found for this gene.

抗原: Synthesized peptide derived from the C-terminal region of human NEIL1

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Ubiquitous.