

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NEIL3 Polyclonal Antibody

多克隆抗体

产品货号: PA05247

第三版

描述:NEIL3 belongs to a class of DNA glycosylases homologous to the bacterial Fpg/Nei family. These glycosylases initiate the first step in base excision repair by cleaving bases damaged by reactive oxygen species and introducing a DNA strand break via the associated lyase reaction. NEIL3 (Nei Like DNA Glycosylase 3) is a Protein Coding gene. Among its related pathways are Telomere C-strand (Lagging Strand) Synthesis. GO annotations related to this gene include nucleic acid binding and single-stranded DNA binding. An important paralog of this gene is NEIL2.

抗原:Synthesized peptide derived from the C-terminal region of human NEIL3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Detected in thymus and testis. Expressed at higher levels in tumor tissue than in the corresponding normal tissues except for pancreas and testis where expression is higher in normal than tumor tissue.