

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NEURL2 Polyclonal Antibody

多克隆抗体

产品货号: PA05259

第三版

描述:This gene encodes a protein that is involved in the regulation of myofibril organization. This protein is likely the adaptor component of the E3 ubiquitin ligase complex in striated muscle, and it regulates the ubiquitin-mediated degradation of beta-catenin during myogenesis. Alternative splicing of this gene results in multiple transcript variants. NEURL2 (Neutralized E3 Ubiquitin Protein Ligase 2) is a Protein Coding gene. An important paralog of this gene is NEURL4.

抗原:Synthesized peptide derived from the Internal region of human Neuralized-2

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed specifically in skeletal and cardiac muscles.