



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat NRSN1 Polyclonal Antibody

多克隆抗体

产品货号: PA05260

第三版

描述:NRSN1 (Neurensin 1) is a Protein Coding gene. Diseases associated with NRSN1 include Louse-Borne Relapsing Fever. An important paralog of this gene is NRSN2.

抗原:Synthesized peptide derived from the Internal region of human Neurensin-1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in brain. Not detectable in other tissues tested.