

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NCAN Polyclonal Antibody

多克隆抗体

产品货号: PA05264

第三版

**描述:** Neurocan is a chondroitin sulfate proteoglycan thought to be involved in the modulation of cell adhesion and migration. May modulate neuronal adhesion and neurite growth during development by binding to neural cell adhesion molecules (NG-CAM and N-CAM). Chondroitin sulfate proteoglycan; binds to hyaluronic acid.

**抗原:** Synthesized peptide derived from the Internal region of human Neurocan

**配方:**

**如何使用:** 加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:** PBS (pH7.4) ， 1% BSA

**应用:** WB 1 ~ 5  $\mu$ g/ml.

**特异性:** Brain.