

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat NF1 Polyclonal Antibody

多克隆抗体

产品货号: PA05265

第三版

描述: This gene product appears to function as a negative regulator of the ras signal transduction pathway. Mutations in this gene have been linked to neurofibromatosis type 1, juvenile myelomonocytic leukemia and Watson syndrome. The mRNA for this gene is subject to RNA editing (CGA>UGA->Arg1306Term) resulting in premature translation termination. Alternatively spliced transcript variants encoding different isoforms have also been described for this gene. NF1 (Neurofibromin 1) is a Protein Coding gene. Diseases associated with NF1 include Neurofibromatosis-Noonan Syndrome and Neurofibromatosis, Type 1. Among its related pathways are Development VEGF signaling via VEGFR2 - generic cascades and RET signaling. GO annotations related to this gene include binding and phosphatidylcholine binding.

抗原: Synthesized peptide derived from the Internal region of human Neurofibromin

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Detected in brain, peripheral nerve, lung, colon and muscle.