

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat CSPG5 Polyclonal Antibody

多克隆抗体

产品货号: PA05266

第三版

**描述:**The protein encoded by this gene is a proteoglycan that may function as a neural growth and differentiation factor. Several transcript variants encoding different isoforms have been found for this gene. CSPG5 (Chondroitin Sulfate Proteoglycan 5) is a Protein Coding gene. Diseases associated with CSPG5 include Spontaneous Ocular Nystagmus and Kabuki Syndrome 1. Among its related pathways are Chondroitin sulfate/dermatan sulfate metabolism and Glycosaminoglycan metabolism. GO annotations related to this gene include growth factor activity.

**抗原:**Synthesized peptide derived from the Internal region of human Neuroglycan C

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$  g/ml.

**特异性:**Restricted to brain (at protein level).