

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NMS Polyclonal Antibody

多克隆抗体

产品货号: PA05267

第三版

描述: This gene encodes a member of the neuromedin family of neuropeptides. The encoded preproprotein is proteolytically processed to generate a biologically active neuropeptide that plays a role in the regulation of circadian rhythm, anorexigenic action, antidiuretic action, cardiovascular function and stimulation of oxytocin and vasopressin release.

抗原: Synthesized peptide derived from the C-terminal region of human Neuromedin-S

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: