

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NMU Polyclonal Antibody

多克隆抗体

产品货号: PA05268

第三版

描述: This gene encodes a member of the neuromedin family of neuropeptides. The encoded protein is a precursor that is proteolytically processed to generate a biologically active neuropeptide that plays a role in pain, stress, immune-mediated inflammatory diseases and feeding regulation. Increased expression of this gene was observed in renal, pancreatic and lung cancers. Alternative splicing results in multiple transcript variants encoding different isoforms. Some of these isoforms may undergo similar processing to generate the mature peptide. NMU (Neuromedin U) is a Protein Coding gene. Diseases associated with NMU include Trachea Leiomyoma and Squamous Cell Papilloma. Among its related pathways are Peptide ligand-binding receptors and RET signaling. GO annotations related to this gene include receptor binding and neuromedin U receptor binding.

抗原: Synthesized peptide derived from the C-terminal region of human Neuromedin-U

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: Expressed throughout the enteric nervous system with highest levels being found in the jejunum.