

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NTM Polyclonal Antibody

多克隆抗体

产品货号: PA05272

第三版

描述: This gene encodes a member of the IgLON (LAMP, OBCAM, Ntm) family of immunoglobulin (Ig) domain-containing glycosylphosphatidylinositol (GPI)-anchored cell adhesion molecules. The encoded protein may promote neurite outgrowth and adhesion via a homophilic mechanism. This gene is closely linked to a related family member, opioid binding protein/cell adhesion molecule-like (OPCML), on chromosome 11. Multiple transcript variants encoding different isoforms have been found for this gene. NTM (Neurotrimin) is a Protein Coding gene. Diseases associated with NTM include Childhood Pilocytic Astrocytoma and Juvenile Pilocytic Astrocytoma. Among its related pathways are Post-translational modification- synthesis of GPI-anchored proteins and Metabolism of proteins. An important paralog of this gene is OPCML.

抗原: Synthesized peptide derived from the C-terminal region of human Neurotrimin

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μ g/ml.

特异性: