

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ANO7 Polyclonal Antibody

多克隆抗体

产品货号: PA05279

第三版

描述: Does not exhibit calcium-activated chloride channel (CaCC) activity. May play a role in cell-cell interactions.

抗原: Synthesized peptide derived from the C-terminal region of human NGEF

配方:

如何使用: 加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液: PBS (pH7.4) ， 1% BSA

应用: WB 1 ~ 5 μg/ml.

特异性: Specifically expressed in epithelial cells of the prostate (at protein level).