

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Monkey NKX6-3 Polyclonal Antibody

多克隆抗体

产品货号: PA05293

第三版

描述:The NKX family of homeodomain proteins controls numerous developmental processes. Members of the NKX6 subfamily, including NKX6-3, are involved in development of the central nervous system (CNS), gastrointestinal tract, and pancreas. NKX6-3 (NK6 Homeobox 3) is a Protein Coding gene. GO annotations related to this gene include sequence-specific DNA binding. An important paralog of this gene is NKX6-2.

抗原:Synthesized peptide derived from the Internal region of human Nkx-6.3

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4)，1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性: