

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NMUR1 Polyclonal Antibody

多克隆抗体

产品货号: PA05296

第三版

**描述:**NMUR1 (Neuromedin U Receptor 1) is a Protein Coding gene. Among its related pathways are Peptide ligand-binding receptors and RET signaling. GO annotations related to this gene include G-protein coupled receptor activity and neuromedin U receptor activity. An important paralog of this gene is NMUR2.

**抗原:**Synthesized peptide derived from the N-terminal region of human NMUR1

**配方:**

**如何使用:**加1ml超纯水重溶

**稳定性:** -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

**稀释液:**PBS (pH7.4) ， 1% BSA

**应用:**WB 1 ~ 5  $\mu$ g/ml.

**特异性:**Expressed in greatest abundance in peripheral organs, particularly in elements of the gastrointestinal and urogenital systems with highest levels in testes. In central nervous system structures express levels are much lower than those seen in peripheral organs. Within the CNS, has been detected in highest abundance in the cerebellum, dorsal root ganglia, hippocampus, and spinal cord.