

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NOL8 Polyclonal Antibody

多克隆抗体

产品货号: PA05299

第三版

描述:NOL8 binds Ras-related GTP-binding proteins (see MIM 608267) and plays a role in cell growth (Sekiguchi et al., 2004 [PubMed 14660641]). Plays an essential role in the survival of diffuse-type gastric cancer cells. Acts as a nucleolar anchoring protein for DDX47. May be involved in regulation of gene expression at the post-transcriptional level or in ribosome biogenesis in cancer cells.

抗原:Synthesized peptide derived from the C-terminal region of human Nop132

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Expressed in various diffuse-type gastric cancers. Detected at lower levels in skeletal muscle.