
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat NPDC1 Polyclonal Antibody

多克隆抗体

产品货号: PA05307

第三版

描述:NPDC1 (Neural Proliferation, Differentiation And Control 1) is a Protein Coding gene.

抗原:Synthesized peptide derived from the Internal region of human NPDC-1

配方:

如何使用:加1ml超纯水重溶

稳定性: -20 ° C保存条件下，冻干粉,保质期为五年；液体，保质期为两年。

稀释液:PBS (pH7.4) ， 1% BSA

应用:WB 1 ~ 5 μ g/ml.

特异性:Strongly expressed in adult brain; especially in hippocampus, frontal lobe and temporal lobe.